

AddGenotype qPCR Master (2x conc.)

Product Information

Product Code 78520

Provided with AddGenotype qPCR Master (2x conc.) 1.0 mL x 5 tubes

Storage Conditions -10°C ~ -30°C

Stability Stable for 2 years from manufacturing date.

Components of AddGenotype qPCR Master (2x conc.)

AddStart Taq DNA Polymerase, Tris-HCl (pH8.8), Potassium Chloride, Ammonium Sulfate, MgCl₂, Protein stabilizer and dNTP mixture.

Description

AddGenotype qPCR Master is a ready - to - use mixture with AddStart Taq DNA Polymerase, reaction buffer with MgCl₂, dNTPs, PCR enhancer and dNTP mixture for TaqMan Probe based real - time PCR.

Especially this Master is optimized qPCR Master for genotyping by TaqMan probe based qPCR.

Applications

- Genotype Analysis by TaqMan Probe based Real - time PCR

Nucleic Acid Amplification Procedure

1. Add the following components to a thin-walled PCR tube:

| | |
|------------------------------------|-------------|
| Nuclease-Free Water | x µl |
| AddGenotype qPCR Master (2x conc.) | 10 µl |
| Forward primer (10 µM) | 0.25 ~ 2 µl |
| Reverse primer (10 µM) | 0.25 ~ 2 µl |
| TaqMan Probe (10 µM) | 0.25 ~ 2 µl |
| (Optional) 50x ROX dye | x µl |
| DNA template | x µl |
| Total reaction volume | 20 µl |

2. PCR cycling

| | |
|-----------------------------------|---|
| Initial denaturation | 95°C, 5 - 10 min |
| PCR cycling (35 - 45 cycles) | 95°C, 10 - 30 sec 60 - 65°C, 30 - 60 sec |

[Note] 50x ROX dye

ROX dye can be included in the reaction to normalize the fluorescent reporter signal, for instruments which are compatible with that option. 50x ROX is a 25 µM concentration. Use the following table to determine the amount of ROX to use with a particular instrument.

| Instrument | Final ROX concentration |
|---|-------------------------|
| AB 7000, 7300, 7700, 7900HT, 7900 Fast, StepOne and StepOnePlus | 500 nM |
| AB 7500, 7500 Fast, Stratagene Mx3000P, Mx3005P and Mx4000 | 50 nM |

Manufacture

에드바이오메디텍

www.addbiomeditek.com / addbiomeditek@addbiomeditek.com