

AddStart Taq DNA Polymerase

Product Information

- Product Code** 18500
- Provided with**
1. AddStart Taq DNA Polymerase (5 U/μl) 100 μl
 2. 10x Reaction Buffer 1.0 mL
 3. 10 mM dNTP Mix. 1.0 mL
- Storage Conditions** -10°C ~ -30°C
- Stability** Stable for 2 years from manufacturing date.
- Storage Buffer** 20 mM Tris-HCl (pH 8.0), 100 mM KCl, 3 mM MgCl₂, 1 mM DTT, 0.1% Tween® 20 and 50% (v/v) glycerol
- 10 mM dNTP Mix** Premixed aqueous solutions of dATP, dCTP, dGTP and dTTP, each at a final concentration of 2.5 mM

Description

AddStart Taq DNA Polymerase is a highly thermostable recombinant DNA polymerase derived from the thermophile, *Thermus aquaticus*, and is a AddStart Taq DNA Polymerase by specific anti - Taq monoclonal antibody.

AddStart Taq DNA polymerase catalyzes the 5' → 3' synthesis of DNA but has no detectable 3' → 5' proofreading exonuclease activity, and possesses low 5' → 3' exonuclease activity, which results in a 3'-dA overhang on the PCR product. Especially, this enzyme can be applied to multiplex PCR, allele specific PCR, SNP analysis and real-time PCR by fluorescent intercalating dye like SYBR Green I® and TaqMan® Probe.

Applications

- Real-time PCR
- General PCR
- Colony PCR
- Screening PCR
- Cloning

Nucleic Acid Amplification Procedure

1. Add the following components to a thin-walled PCR tube:

Nuclease-Free Water	x μl
10x Reaction Buffer	2 μl
10mM dNTP Mixture	2 μl
Forward primer (10 μM)	0.25 ~ 2 μl
Reverse primer (10 μM)	0.25 ~ 2 μl
DNA template	x μl
AddStart Taq DNA Polymerase (5 U/μl)	0.2 μl
Total reaction volume	20 μl

* Recommendation of template DNA concentration in a 20 μl reaction volume

- 1) Human genomic DNA: 0.1 ng ~ 1 μg
- 2) Bacterial, Plant genomic DNA: 0.1 ng ~ 100 ng
- 3) Plasmid DNA: 0.01 ng ~ 5 ng

2. PCR cycling

Initial denaturation	95°C, 5 - 10 min
PCR cycling (25 - 40 cycles)	95°C, 15 - 30 sec
	55 - 65°C, 15 - 30 sec
	72°C, 30 - 60 sec
Final extension	72°C, 5 min
Hold	12°C, ∞

Manufacture

애드바이오메디텍

www.addbiomeditek.com / addbiomeditek@addbiomeditek.com