

# Add Taq DNA Polymerase

## Product Information

<b>Product Code</b>	17500
<b>Provided with</b>	1. Add Taq DNA Polymerase (5 U/μl) 100 μl 2. 10x Reaction Buffer 1.0 mL 3. 10 mM dNTP Mix. 1.0 mL
<b>Storage Conditions</b>	-10°C ~ -30°C
<b>Stability</b>	Stable for 2 years from manufacturing date.
<b>Storage Buffer</b>	20 mM Tris-HCl (pH 8.0), 100 mM KCl, 3 mM MgCl <sub>2</sub> , 1 mM DTT, 0.1% Tween® 20 and 50% (v/v) glycerol
<b>10 mM dNTP Mix</b>	Premixed aqueous solutions of dATP, dCTP, dGTP and dTTP, each at a final concentration of 2.5 mM

## Description

Add Taq DNA Polymerase is a highly thermostable recombinant DNA polymerase derived from the thermophile, *Thermus aquaticus*. The molecular weight of the recombinant protein is 94kD. Add Taq DNA polymerase catalyzes the 5' → 3' synthesis of DNA but has no detectable 3' → 5' proofreading exonuclease activity, and possesses low 5' → 3' exonuclease activity, which results in a 3'-dA overhang on the PCR product. Normal maximum amplicon size is ~5kb and the polymerization rate is 1 kb per 1 min., but maximum rate is 2 kb per 1 min.

## Applications

- General PCR
- Colony PCR
- Screening PCR
- Cloning

## Nucleic Acid Amplification Procedure

1. Add the following components to a thin-walled PCR tube:

Nuclease - Free Water	x μl
10x Reaction Buffer	2 μl
10mM dNTP Mixture	2 μl
Forward primer (10 μM)	0.25 ~ 2 μl
Reverse primer (10 μM)	0.25 ~ 2 μl
DNA template	x μl
Add Taq DNA Polymerase (5 U/μl)	0.2 μl
Total reaction volume	20 μl

\* Recommendation of template DNA concentration in a 20 μl reaction volume

- 1) Human genomic DNA: 0.1 ng ~ 1 μg
- 2) Bacterial, Plant genomic DNA: 0.1 ng ~ 100 ng
- 3) Plasmid DNA: 0.01 ng ~ 5 ng

2. PCR cycling

Initial denaturation	95°C, 5 min
PCR cycling ( 25 - 40 cycles )	95°C, 15 - 30 sec 55 - 65°C, 15 - 30 sec 72°C, 30 - 60 sec
Final extension	72°C, 5 min
Hold	12°C, ∞

## Manufacture

애드바이오메디텍

www.addbiomeditek.com / addbiomeditek@addbiomeditek.com