

AddProbe RT-qPCR Kit and Master with UDG

Product Information

Product Code 74201, 74201U

Provided with

Cat. No. 74201: AddProbe RT-qPCR Kit with UDG 1. 20x Enzyme Solution 0.1 mL
2. 5x Reaction Buffer 0.4 mL

Cat. No. 74201U: AddProbe RT-qPCR Master with UDG (2x conc.) 1.0 mL

Storage Conditions -10°C ~ -30°C

Stability Stable for 1 year from manufacturing date.

Components of AddProbe RT-qPCR Kit and Master with UDG

AddStart Taq DNA Polymerase, Reverse Transcriptase, RNase Inhibitor, UDG, Tris - HCl (pH8.5), Potassium Chloride, Ammonium Sulfate, MgCl₂, DTT, Protein stabilizer and dNTP mixture.

Description

AddProbe RT-qPCR Kit and Master with UDG provides sensitive and easy-to-use components which contain all the reagents for first strand cDNA synthesis and PCR reaction in one-tube by using TaqMan® Probe, and is designed for high sensitivity and specificity on various real-time instruments. Especially, thermostable MMLV RTase (RNase H-), AddStart Taq DNA Polymerase, RNase Inhibitor and heat-labile UDG are included in 20x Enzyme Solution.

AddProbe RT-qPCR Master with UDG (2x conc.) is a ready-to-use mixture included all components for one-step RT-qPCR in one tube.

Applications

- TaqMan Probe based Real-time PCR for RNA template
- Quantitative real-time PCR
- Infectious disease detection
- RNA virus detection

Manufacture

애드바이오메디텍

www.addbiomeditek.com / addbiomeditek@addbiomeditek.com



ADDBIOMEDITEK Co., Ltd.

Nucleic Acid Amplification Procedure

1-1. Add the following components to a thin-walled PCR tube for AddProbe RT-qPCR Kit:

Nuclease-Free Water	x µl
5x Reaction Buffer	4.0 µl
Forward primer (10 µM)	0.25 ~ 2.0 µl
Reverse primer (10 µM)	0.25 ~ 2.0 µl
TaqMan Probe (10 µM)	0.25 ~ 2.0 µl
(Optional) 50x ROX dye	x µl
RNA template	x µl
20x Enzyme Solution	1.0 µl
Total reaction volume	20 µl

1-2. Add the following components to a thin-walled PCR tube for AddProbe RT-qPCR Master:

Nuclease-Free Water	x µl
2x AddProbe RT-qPCR Master (2x conc.)	10.0 µl
Forward primer (10 µM)	0.25 ~ 2.0 µl
Reverse primer (10 µM)	0.25 ~ 2.0 µl
TaqMan Probe (10 µM)	0.25 ~ 2.0 µl
(Optional) 50x ROX dye	x µl
RNA template	x µl
Total reaction volume	20 µl

2. PCR cycling

cDNA synthesis	50°C, 20 min
Initial denaturation	95°C, 5 - 10 min
PCR cycling (35 - 45 cycles)	95°C, 10 - 30 sec 60 - 65°C, 30 - 60 sec

[Note] 50x ROX dye

ROX dye can be included in the reaction to normalize the fluorescent reporter signal, for instruments which are compatible with that option. 50x ROX is a 25 µM concentration. Use the following table to determine the amount of ROX to use with a particular instrument.

Instrument	Final ROX concentration
AB 7000, 7300, 7700, 7900HT, 7900 Fast, StepOne and StepOnePlus	500 nM
AB 7500, 7500 Fast, Stratagene Mx3000P, Mx3005P and Mx4000	50 nM